GB

Frog Appliance*

Laboratory instructions

D

Frosch Apparatur*

Laborinstruktionen

F

L'appareil Grenouille*

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Phone +49 (0) 7231 459-0 · Fax +49 (0) 7231 459-102

info@forestadent.com · www.forestadent.com

FORESTADENT[®]
GERMAN PRECISION IN ORTHODONTICS

1

**Preliminary laboratory requirements / Laborbedarf /
Etape préliminaire en laboratoire**

GB

- A) Maxillary molar bands with lingual sheaths cemented in place.
- B) An accurate model of the maxillary arch including the palatal tissue, the occlusal and lingual surfaces of the teeth, and the lingual sheaths of the first molar bands.

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2a

GB

- Apply a separating medium to the model.
- Bend anchor wires from 0,7mm/.028" FORESTANIT. The anchor wires should lie in the embrasures distal to the anchor teeth.
- The FROG screw is placed on the model with the distal of the appliance being aligned antero-posteriorly with the mesial of the lingual sheaths.

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
3

- GB** d) Occluso-apically the FROG should be placed **approximately 10 mm to 12 mm from the occlusal surface**. This will place the appliance at approximately the center of resistance of the molars for bodily tooth movement. **This is extremely important!** If placed too far occlusally, distal crown tipping will occur. If placed too far apically the roots will tip distally.
- e) Make sure the angle of the screw driver to the FROG screw does not exceed 15° when placed inside the FROG screw for activation.

Caution! Only use the original screw driver that is delivered with the FROG!

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
4

GB

Caution: Isolate the FROG-screw before applying acrylic and make sure there is no acrylic flowing inside the screw mechanism!

- f) Fabricate a Nance button with the anterior extensions of the FROG embedded in it.
- g) The body of the FROG can be removed from the button to allow easier polishing.

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**GB**

h) Use the preformed .032" stainless steel blank to customize the distalizing spring (FROGLEGS).

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6

GB

i) The polished appliance is secured with an elastic and is ready for delivery to the patient.

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