# Track<sup>®</sup> PC

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**Chemical characterization:**

PC: Polycarbonate. Non-hazardous to health, BfR- listed (German Federal Institute for Risk Assessment) thermoplastic foil with certified biocompatibility, insoluble in water, inactive, no risk to ground water.

**Disposal/Recycling:**

Domestic waste (reduced, air pollutant incineration) / unmixed plastic recycling.

**Area of application:**

Essix retainers, miniplast splints, bite guards, retention splints, adapters.

**Material properties:**

Tough, clear foils with optimal transparency, fracture resistant, abrasion resistant and exceptionally impact resistant. Excellent vacuum-forming properties using the press technique, polycarbonate bonds with cold-curing resin after activation with a bonding agent. Foils should be dried before vacuum forming to avoid bubbling. Dry the foils in a drying chamber at 80°C – 100°C for approx. 4 hours. Dried foils can be vacuum formed for several days. Foils can be dried again without affecting the quality.

**Preparation:**

Polish to a high-lustre polish with acrylic polishers. Apply a little monomer beforehand to sections that are to be built up with acrylic.

**Plastification:**

(For units without heating time control) Place the foil in the unit, test the softness with an instrument, then vacuum form, ensure the plastification temperature is not too high (approx. 180°C).

**Pack contents:** Track® PC, transparent, Ø 125 mm

Description	Thickness	Content	Best.-Nr.	Heating time for MINISTAR® / Biostar® by SCHEU-DENTAL	
				Heating time	Code
Track® PC	1,0 mm	100	408-0700	ca. 50 sec	164

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