

## Die Frosch-Apparatur – zur einfachen Molarendistalisation\* Frog Appliance – the Simplified Molar Distalizer\*

Info  
N° 177



### Die Frosch-Apparatur






Bei Klasse-II-Malokklusionen – und manchmal auch bei Non-Extraktionsbehandlungen von OK-Engständen – ist das Ziel oftmals die Distalisation der OK-Molaren.

Die Frosch-Apparatur bietet eine hervorragende Molarendistalisation und körperliche Zahnbewegung ohne Kippungen. Einbau und Aktivierung der Frosch-Apparatur sind ebenso einfach wie das Herausnehmen und die Nachjustierung der Feder. Die Apparatur ist von der Frontseite her leicht zu aktivieren, indem man die Schraube entgegen dem Uhrzeigersinn dreht.

### The Frog Appliance

Distalizing maxillary first molars is often an objective in treatment plans involving Class II malocclusions and it is sometimes indicated for non-extraction treatments with maxillary crowding.

The Frog appliance offers outstanding molar distalization without crown or root tipping. The assembling and activation of the Frog is easy. Even the spring is easily removed and adjusted. The appliance is easily activated from the anterior by simply turning the screw counter clockwise.

Order No. CE 0297		Beschreibung description
169-1000		<b>Frosch-Apparatur – Kit für 1 Patienten</b> 1 Schraube, 1 vorgeformte Feder, 1 Schraubendreher, 1 Labor- und Einbauanleitung <b>FROG-appliance – kit for 1 patient</b> 1 screw, 1 preformed spring, 1 screwdriver 1 Laboratory and delivery instruction
169-5000		<b>Frosch-Apparatur – Kit für 5 Patienten</b> ohne Schraubendreher 5 Schrauben, 5 vorgeformte Federn, 1 Labor- und Einbauanleitung <b>FROG-appliance – kit for 5 patients</b> without screwdriver 5 screws, 5 preformed springs, 1 laboratory and delivery instruction
169-1523		<b>Schraube für Frosch Apparatur</b> (1 Pack = 5 Stück) <b>Screw for FROG-appliance</b> (1 pack = 5 pcs.)
307-1523		<b>Vorgeformte Feder für Frosch Apparatur</b> (1 Pack = 5 Stück) <b>Preformed spring for FROG-appliance</b> (1 pack = 5 pcs.)
100-1523 CE		<b>Schraubendreher für Frosch-Apparatur</b> (1 Pack = 5 Stück) <b>Screwdrivers for FROG-appliance</b> (1 pack = 5 pcs.)