




Ni-Free Products

Bukkalröhrchen zum Kleben – Nickelfrei Buccal tubes for direct bonding – Ni-Free

CE 0297



Pad sandgestrahlt
Sandblasted pad

OK / Maxillary	Zahn tooth	Torque	Angulation	distaler Offset distal offset	auf/on Pad (10 Stück pro Packung/pcs. per pack)		
					rechts/right	links/left	
occlusal headgear tube  convertible	$\frac{6}{ } \frac{6}{ }$	-14°	0°	14°	.018" x .025"	741M0704	740M0704
					.022" x .028"	741M0702	740M0702
	$\frac{7}{ } \frac{7}{ }$	-14°	0°	14°	.018" x .025"	741M0804	740M0804
					.022" x .028"	741M0802	740M0802
UK / Mandibular	$\frac{6}{ } \frac{6}{ }$	-30°	0°	4°	.018" x .025"	741M1704	740M1704
					.022" x .028"	741M1702	740M1702
	$\frac{7}{ } \frac{7}{ }$	-30°	0°	4°	.018" x .025"	741M1804	740M1804
					.022" x .028"	741M1802	740M1802